



POLYSULFONE DIALYZER



# RETHINKING HEMODIALYSIS



2025

[www.NephroCan.com](http://www.NephroCan.com)



## CANADIAN LEADER IN HEMODIALYSIS



As a trusted Canadian provider, we are committed to delivering high-quality hemodialysis consumables, machinery, and equipment worldwide. Our goal is to ensure that our offerings are accessible and enhance patient care. We are dedicated to pioneering innovative therapies that revolutionize renal therapy, improving outcomes and enhancing the quality of life for patients globally.

# HIGH FLUX POLYSULFONE HEMODIALYZER FILTERS

Polysulfone High Flux Hollow Fiber Hemodialyzer Specifications*								
Sterilized with Gamma Irradiation (R)								
	Blood Flow Rate (mL/min)	Nephro PS100R	Nephro PS130R	Nephro PS140R	Nephro PS160R	Nephro PS180R	Nephro PS200R	Nephro PS220R
Clearances (mL/min)	Urea	QB= 200	190	195	196	197	198	199
		QB= 300	258	264	272	284	286	293
	Creatinine	QB= 200	173	180	186	193	194	198
		QB= 300	217	234	247	267	276	291
	Phosphate	QB= 200	160	174	179	188	192	196
		QB= 300	199	219	230	249	259	288
	Vitamin B <sub>12</sub>	QB= 200	113	130	139	153	161	168
		QB= 300	127	146	159	178	190	221
	Inulin	QB= 200	78	93	104	117	125	135
		QB= 300	84	103	113	130	141	163
Surface Area (m <sup>2</sup> )		1.0	1.3	1.4	1.6	1.8	2.0	2.2
Blood Priming Volume (mL)		59	69	74	86	105	107	116
Ultrafiltration Coefficient (mL/mmHg.hr)		34.6	46.4	50.6	59.4	63.7	74.6	80.9
KoA (mL/min)		930	1030	1190	1570	1660	1900	2160
Membrane								
Membrane Material		Dry Non-Cellulosic Membrane, High Flux Polysulfone						
Wall Thickness (μm)		40 ± 7						
Inner Diameter (μm)		200 ± 15						
Maximum TMP (mmHg)		500						
Sieving Coefficient**								
β2-Microglobulin		> 0.8						
Albumin		< 0.001						

QB: Blood flow rate


QF: Ultrafiltration rate

QD: Dialysate flow rate

\*Specifications and performance data at QB = 200/300 mL/min, QD = 500 mL/min, QF = 0 mL/min, T: 37 °C. Performance data was measured in vitro according to standards BS EN ISO 8637-1:2020. Clearance data may vary depending on testing conditions

\*\*Sieving coefficient is an absolute value that is related to the hollow membrane fiber

# REMOVAL OF UREMIC TOXINS IN HIGH FLUX DIALYZERS



## Dialyzer Sizes, High Flux, Gamma Sterilization

Specifications and performance based on NephroCan's PS 1.8M<sup>2</sup> High Flux Polysulfone membranes

QB = 300 mL/min, QD = 500 mL/min, QF = 0 mL/min

\*Specifications and performance based on public data available on average competing hemodialyzers with Polysulfone membranes  
QB = 300 mL/min, QD = 500 mL/min, QF = 0 mL/min

## NEPHROFILTER CHARACTERISTICS



# LOW FLUX POLYSULFONE HEMODIALYZER FILTERS

# REMOVAL OF UREMIC TOXINS IN LOW FLUX DIALYZERS


Polysulfone Low Flux Hollow Fiber Hemodialyzer Specifications*									
Sterilized with Gamma Irradiation (R)									
		Blood Flow Rate (mL/min)	Nephro PS10R	Nephro PS13R	Nephro PS14R	Nephro PS16R	Nephro PS18R	Nephro PS20R	Nephro PS22R**
Clearances (mL/min)	Urea	QB= 200	187	195	196	197	198	199	199
		QB= 300	238	250	253	274	284	285	291
	Creatinine	QB= 200	167	180	184	189	189	191	193
		QB= 300	201	223	226	243	265	268	283
	Phosphate	QB= 200	143	153	157	164	174	177	182
		QB= 300	161	182	188	196	215	220	234
	Vitamin B <sub>12</sub>	QB= 200	79	93	98	120	139	145	149
		QB= 300	84	100	117	135	156	161	176
	Inulin	QB= 200	-	-	-	-	-	-	-
		QB= 300	-	-	-	-	-	-	-
Surface Area (m <sup>2</sup> )		1.0	1.3	1.4	1.6	1.8	2.0	2.2	
Blood Priming Volume (mL)		59	69	74	86	105	107	110	
Ultrafiltration Coefficient (mL/mmHg.hr)		7.3	9.5	11.9	13.9	16.7	17.8	19.8	
KoA (mL/min)		698	824	977	1240	1570	1610	1980	
Membrane									
Membrane Material		Dry Non-Cellulosic Membrane Low Flux Polysulfone							
Wall Thickness (μm)		40 ± 7							
Inner Diameter (μm)		200 ± 15							
Maximum TMP (mmHg)		500							

QB: Blood flow rate

QF: Ultrafiltration rate

QD: Dialysate flow rate

\*Specifications and performance data at QB = 200/300 mL/min, QD = 500 mL/min, QF = 0 mL/min, T: 37 °C. Performance data was measured in vitro according to standards BS EN ISO 8637-1:2020. Clearance data may vary depending on testing conditions



## Dialyzer Sizes, Low Flux, Gamma Sterilization

Specifications and performance based on NephroCan's PS 1.8M<sup>2</sup> Low Flux Polysulfone membranes

QB = 300 mL/min, QD = 500 mL/min, QF = 0 mL/min

\*Specifications and performance based on public data available on average competing hemodialyzers with Polysulfone membranes

QB = 300 mL/min, QD = 500 mL/min, QF = 0 mL/min

## NEPHROFILTER CHARACTERISTICS





**NephroCan Inc.**

Suite 940-1040 West Georgia St. Vancouver, BC V6E 4H1  
Canada



**NephroCan S.r.l.**

Via Savoia 78 00198 Roma, Italy  
[www.NephroCan.com/EU-Rep](http://www.NephroCan.com/EU-Rep)



**Phone:**

+1 (604) 922-3110

**Fax:**

+1 (604) 922-3160

**Email:**

[info@NephroCan.com](mailto:info@NephroCan.com)  
[sales@NephroCan.com](mailto:sales@NephroCan.com)

**Website:**

[www.NephroCan.com](http://www.NephroCan.com)

