NephroCan

NephroChair Characteristics	
Max. patient weight	200 kg
Emergency position	Trendelenburg position in handset and foot pedal Flat position in foot pedal
Number of actuators	Fully adjustable (back rest, seat rest, leg rest, foot rest, and height) and equipped with emergency settings that function in all chair positions
Back rest angle	0° to 105°
Foot rest	0° to 57°
Seat rest	0° to 12°
Input voltage of control box	220V-230 V 50 Hz
Output voltage of control box	24 V DC
Input voltage of handset, actuator and foot pedal	24 V DC
Output voltage of handset, actuator and foot pedal	24 V DC
Main voltage	220V-230 V 50 Hz
Power rating	0.7 kW
Electrical protection	Class I
Degree of protection against penetration of solids / liquids	IPX4
Work environment and storage condition	5 to 45 °C
Mattress and cover	Memory foam mattress, 90°C washable polyurethane (PU) cover
Standards	EN 60601-1, EN 60601-1-2, EN ISO 13485, 93/42/EEC Directive
Options	Reading list, foot pedal and colour palette
OTHER FEATURES	
Mobile, lightweight, and compact with rounded corners for safety	
Equipped with patient weighing scale	
Highly customizable to fit any design requirements for centers globally	
Hand control unit and central lockable casters	
ABS thermoplastic frame cover providing rigidity and impact (heat resistance)	
Compliant with all safety test requirements, including electromedical standards	

No representation or warranty, expressed or implied, is given as to the correctness or completeness of the information provied in this document, and NephroCan assumes no liability. The information and content provided is private intellectual property of NephroCan and only to be used for purposes as intended by NephroCan. The content of this document shall not be altered, changed, or replicated. NephroCan claims possessory rights over the content of this document.